

Restore


Enhancing Comfort and Care
in Rehabilitation Ensuites

Custodian® Shower Rose
 Safe and secure design providing a comfortable showering experience. Unique design allows for easy room adaptation to accommodating Carer's hand showers with no additional diverter required.

The Restore Ensuite

Enhancing comfort, safety and dignity for users and carers while enabling independence through intelligent design.

Removable Shower Head



RESTORE TOILET FLUSH VALVE



RH413S

RESTORE SHOWER CONTROL PANEL



RHMIXBT

RESTORE BASIN CONTROL PANEL



RHBSNMIXBT

SHOWER ROSE & REMOVABLE SHOWER HEAD



CNHNDSHHS-D

CARE 800 BACK-TO-WALL TOILET PAN




ENCARE800BTW

KOLO NOVA PRO MENTAL HEALTH KIT



KNP311510-100MHKIT



Example of outlet self draining pipework


Stagnation Management, Self Draining and Programmable Flushing

The restore intelligent control hub is the ideal system for stagnation management as regular pipework flushing can be enabled for every fixture.

The Restore smartboard can be installed above the outlets creating a self draining system further reducing bacterial growth risks in water systems.




To Basin



To Shower

To Toilet



SMARTFLOW® WATER MANAGEMENT

Parameters are preset set within the eTMV Control Hub or are adjustable when connected to a Smart Flow System.

STAGNATION MANAGEMENT

Maximum protection from bacteria and legionella growth requires daily self flushing where usage patterns may indicate risk of bacterial growth.



SMARTBOARD

Controlled by the Restore Smartboard- eTMV, Smart Hub and solenoids are all housed securely in a lockable, tamper resistant cabinet.

Ability to pre-set maximum temperature for special applications (from ambient to 45°C)

WMS-RH-SB001R
(LID SUPPLIED SEPERATE)

TOILET FLUSHING

The secure electronic touch activation panel for full and half flush directly covers the flushing solenoid mounted behind.

SNAKE EYE SECURING SCREWS

Security screws (snake eye) resist tampering or removal without a tool, while allowing simple front-end install and maintenance access.

INTUITIVE ADJUSTMENT

The low resistance temperature handle provides tactile and precise temperature adjustment while the 'touch on' fast start pre-set temperature delivery is ideal for those with limited dexterity or physical disabilities.

Temperature feedback is provided by the backlit temperature handle, changing from blue through yellow to red as the temperature increases.

Secure Clinical Hand Wash

Enhancing safe hand hygiene in mental health environments

ESQX - POINT OF USE
SENSOR OPERATED TMV



ATM621 + ATWM-CNBSN

SENSOR OP TEPID WATER
BASIN OUTLET



DET5071-660

KOLO NOVA PRO
MENTAL HEALTH KIT



KNP311510-100MHKIT

Restore

Safety & Care in Rehabilitation Ensuites

